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1. Introduction

The **DET40/DET75** detectors are used to detect individual particles like electrons, ions and photons and to determine the arrival time of the particle with respect to an external trigger (Time-of-Flight). The effective detection diameters are $> 40\text{mm}$ or $> 75\text{mm}$, respectively. If you have received a different-size MCP please refer to a separate manual for the mounting procedure. Most of the information given will also be valid unless otherwise stated. The **DET40/DET75** product assembly contains an MCP stack and a timing anode. The MCP stack is a chevron configuration (two MCPs) mounted between two ceramic rings. The timing anode consists of a solid metal anode (default) or a coated ceramic substrate. The timing signal can be picked up from the MCP contacts (front or back) or from the timing anode. The so called “transparent timing anode” with a Germanium coating on a ceramic substrate allows for a pick-up of the image charge on a structured electrode pattern (e.g. a Wedge-and-Strip or LC-delay-line anode) behind the detector (available only for **DET40**). It is recommended to pick up the timing from the MCP in this case.

The pick-up of the timing signal (from anode or MCP contact) can be performed according to Figure 1.1. It is recommended to use short cables between the detector and the amplifier or 50 Ohm coaxial cables with proper termination.


Figure 1.1: RC-coupler for getting a fast timing signal from the MCP (here: Z-stack)

The entire detector is bakeable up to 150°C (250°C on request), the metal contacts on the ceramic rings/substrate are suitable for soldering. The operation requires three DC voltages for MCP front and back contacts and for the anode layer and adequate vacuum feedthroughs, signal decouplers and timing electronics. Optional **RoentDek** products for **DET40/DET75** are

- 4-fold MHV feedthrough with signal pickup (**FT4TP**) and mounting flange (**FT4TP/100/150**)
- $2\times 4\text{kV}$ high voltage bias supply (**HV2/4**) and optional **HVZ(-T)** voltage divider.
- Timing amplifier (e.g. **FAMP1+** amplifier) and **CFD1c/1x** or the **ATR19-2** amplifier & CFD unit.

1.1 Characteristics

1.1.1 Physical Characteristics of Detectors DET40/DET75

Mounting Diameter:	80mm/120mm
Mounting height	15mm/20mm
Baking Temperature:	150°C Maximum
Mounting Flange Size: (only with FT4TP100/150)	CF 100/CF150
Height above Flange:	$\approx 100\text{mm}$ (adjustable)
Mounting Diameter:	95mm/146mm

1.1.2 Physical Characteristics of MCPs

# of MCPs in Chevron stack	2
Diameter:	50mm/87mm *
Active Diameter:	42mm/77mm
Quality Diameter:	40mm/75mm
Aspect Ratio (L/D):	80:1
Thickness:	1mm
Pore diameter:	$13\mu\text{m}$
Bias Angle	20°
Open Area Ratio:	70%

* Version **DET40B** or **DET80**: L/D 60:1, $25\mu\text{m}$ pore (1.5mm thick), bias angle 8° , detection grade, OAR $< 60\%$

Operating Temperature Range: -50 to 70°C
Operating Pressure: $< 2 \times 10^{-6}$ mbar

1.1.3 *Electrical Characteristics of Detector*

Electron Gain @ 2400 Volts: 1×10^7 Minimum
Operating Voltage: 2600 Volts Maximum

1.2 List of assembly parts

- 2 ceramic rings, partially metal coated
- timing anode: solid metal anode (default) or coated ceramic substrate.
- 2 micro-channel plates, matched in resistance
- 4 metal spring clamps
- M3 plastic screws with nuts (for pre-assembly of the MCP stack)
- M2 screws and special nuts

1.3 Detector assembly

All parts, especially the MCP and anode should be handled with great care since the surfaces are very sensitive and should never be touched or scratched.

The ceramic rings should not be exposed to exceeding mechanical and thermal stresses and the assembly should take place under clean and dry conditions.

1.3.1 *Preparation:*

- A mesh can be spot-welded or soldered directly onto the front side of the front ceramic ring. The mesh is then placed 1.5mm in front of the MCP surface.
- Prepare the connection cables for the MCP detector and the anode. You need 3 cables for the MCP front, MCP back and anode layer (a 4th if a mesh is used). The cables can be soldered directly onto the metallization of the ceramic rings and the anode. If so, use preferably the metallization strips which are not located on a hole in the ceramic parts. Please take care that no solder or flux is sprayed over the anode structure or the MCPs.
- Otherwise the cables on the anodes can be fixed by screws at the metallization rings located at holes in the ceramic parts. Before mounting be aware which holes should be used for that purpose. For the anode layer one of the screws which are used for fixing the anode can be used. **RoentDek** provides special contact clamps for the MCP (UHV assembly only).
- Clean all parts **except the MCPs and the anode** in an ultrasonic bath.

Warning:

Although it is possible to clean also the anode or the MCP in an ultrasonic bath (isopropanol) for very short time, there is a considerable risk for damage. If cleaning is necessary contact a **RoentDek** agent before.

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


- Aluminium-, copper-, brass- and non-ferrous metal foils and strips as coils and sheets
- Slitting lines for metal foils and strips
- anoxal® anodized oxidized aluminium
- anoxal®E the conducting material
- anoxal®P the physical surface
- anoxal®C the colour surface
- Coils made of aluminium-, copper- and anoxal®E for electrical applications
- Transformers and inductances and supplies

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




KPM 12-50

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
PD-1161

57/60mm | NEMA23/24
 Stepper Motor with
 Controller / Driver
 sensOstep™ Encoder
 Serial Interface

INFO The PANdrive PD57-1161 and PD60-1161 family is a mechatronic solution including a 57mm or 60mm flange motor, a controller board and a sensOstep™ encoder. It can be controlled via serial interface or operated in standalone mode. Power supply, external encoder, interface and I/Os can be connected with JST connectors.

With the advanced stallGuard2™ the motor load can be detected with high resolution. The outstanding coolStep™ technology for sensorless load dependent current control allows efficient motor operation.

The PC based software development environment TMCL-IDE for the Trinamic Motion Control Language (TMCL™) can be downloaded free of charge from the TRINAMIC website.

MAIN CHARACTERISTICS

- ELECTRICAL DATA**
 - nom. 24V DC supply voltage (10V to 30V)
- MOTOR DATA**
 - flange size 57/60mm | NEMA23/24
- INTERFACE**
 - USB, RS232, RS485
 - step&direction interface
 - inputs for ref. & stop switches
 - general purpose I/Os
- FEATURES**
 - up to 256 times microstepping
 - memory for 2048 TMCL™ commands
 - stallGuard2™ sensorless load detection
 - coolStep™ sensorless load dependent current control
 - microPlyer™ 16 to 256 times microstepping interpolation
 - integrated absolute sensOstep™ encoder with 1024 ppr.
 - automatic ramp generation in hardware
 - on the fly alteration of motion parameters
- SOFTWARE**
 - standalone operation using TMCL or remote controlled operation
 - PC-based (Windows) application development software TMCL-IDE downloadable
- OTHER**
 - RoHS compliant

ORDER CODE	DESCRIPTION
PD57-1-1161	0.55 Nm / QMot motor QSH5718-41-28-055
PD57-2-1161	1.01 Nm / QMot motor QSH5718-51-28-101
PD60-3-1161	2.10 Nm / QMot motor QSH6018-65-28-210
PD60-4-1161	3.10 Nm / QMot motor QSH6018-86-28-310
TMCM-1161-CABLE	Cable loom including all necessary cables (single ended)

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